



Product Datasheet

GAS6 antibody (orb125761)


Description	Rabbit Polyclonal to GAS6.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q14393
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from 293 cells.



Western blot analysis of extracts from 293 cells.



Immunohistochemical analysis of formalin-fixed paraffin-embedded (FFPE) tissue sections.