


## Product Datasheet

### GCNT3 antibody (orb125767)

|                            |                                                        |
|----------------------------|--------------------------------------------------------|
| <b>Description</b>         | Rabbit Polyclonal to GCNT3.                            |
| <b>Species/Host</b>        | Rabbit                                                 |
| <b>Reactivity</b>          | Human                                                  |
| <b>Conjugation</b>         | Unconjugated                                           |
| <b>Tested Applications</b> | IF, IHC                                                |
| <b>Immunogen</b>           | Internal                                               |
| <b>Concentration</b>       | 1mg/ml                                                 |
| <b>Note</b>                | For research use only                                  |
| <b>Clonality</b>           | Polyclonal                                             |
| <b>Uniprot ID</b>          | <b>095395</b>                                          |
| <b>Dilution Range</b>      | IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:40000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                        |



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...