



Product Datasheet

GIMAP2 antibody (orb125776)

Description	Rabbit Polyclonal to GIMAP2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9UG22
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from H...



Confocal immunofluorescence analysis of ...