


Product Datasheet

GLB1L3 antibody (orb125784)

Description	Rabbit Polyclonal to GLB1L3.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8NCI6
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...