



## Product Datasheet

### NTF3 Antibody (orb1257921)


|                            |                                                    |
|----------------------------|----------------------------------------------------|
| <b>Description</b>         | NTF3 Antibody                                      |
| <b>Species/Host</b>        | Rabbit                                             |
| <b>Reactivity</b>          | Human, Mouse                                       |
| <b>Conjugation</b>         | Unconjugated                                       |
| <b>Tested Applications</b> | IHC, WB                                            |
| <b>Immunogen</b>           | Recombinant protein of human NTF3                  |
| <b>Target</b>              | NTF3                                               |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3. |
| <b>Form/Appearance</b>     | Liquid                                             |
| <b>Concentration</b>       | batch dependent                                    |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.        |
| <b>Note</b>                | For research use only                              |
| <b>Application notes</b>   | WB: 1:500 - 1:1000IHC: 1:100 - 1:200               |
| <b>Isotype</b>             | IgG                                                |
| <b>Clonality</b>           | Polyclonal                                         |
| <b>MW</b>                  | Observed: 29kDa                                    |
| <b>Uniprot ID</b>          | <b>P20783</b>                                      |
| <b>Dilution Range</b>      | WB: 1:500 - 1:1000IHC: 1:100 - 1:200               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                    |



Western blot analysis of extracts of mou...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...