


Product Datasheet

IL18R1 Antibody (orb1258079)

Description	IL18R1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-320 of human IL18R1 (NP_003846.1).
Target	IL18R1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:1000IF: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 75kDa
Uniprot ID	Q13478
Dilution Range	WB: 1:500 - 1:1000IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of NCI...