



Product Datasheet

EGLN2 Antibody (orb1258208)

Description	EGLN2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human EGLN2 (NP_444274.1).
Target	EGLN2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IP: 1:20 - 1:50
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 43kDa
Uniprot ID	Q96KS0
Dilution Range	WB: 1:500 - 1:2000IP: 1:20 - 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunoprecipitation analysis of 150 ug e...