

PSMC2 Antibody (orb1258351)

Western blot analysis showing PSMC2 expression in Hsp92, SW620, BT-174, OVCAR3, 293T, NCI-H460, Mouse spleen, and Mouse brain. Molecular weight markers are indicated on the left (130kDa, 100kDa, 70kDa, 55kDa, 40kDa, 35kDa, 25kDa). A band for PSMC2 is visible at approximately 50kDa in all lanes except the mouse brain lane.

This micrograph shows a section of skeletal muscle tissue with a normal architecture. The muscle fibers are arranged in a regular, parallel pattern, and the nuclei are clearly visible along the periphery of the fibers.

Immunohistochemistry
of paraffin-embedde...