



## Product Datasheet

### RAB11A Antibody (orb1258790)

|                            |                                                     |
|----------------------------|-----------------------------------------------------|
| <b>Description</b>         | RAB11A Antibody                                     |
| <b>Species/Host</b>        | Rabbit                                              |
| <b>Reactivity</b>          | Human, Mouse, Rat                                   |
| <b>Conjugation</b>         | Unconjugated                                        |
| <b>Tested Applications</b> | IF, IHC, WB                                         |
| <b>Immunogen</b>           | Recombinant protein of human RAB11A.                |
| <b>Target</b>              | RAB11A                                              |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.  |
| <b>Form/Appearance</b>     | Liquid                                              |
| <b>Concentration</b>       | batch dependent                                     |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.         |
| <b>Note</b>                | For research use only                               |
| <b>Application notes</b>   | WB: 1:500 - 1:2000IHC: 1:50 - 1:100IF: 1:50 - 1:100 |
| <b>Isotype</b>             | IgG                                                 |
| <b>Clonality</b>           | Polyclonal                                          |
| <b>MW</b>                  | Observed: 24kDa                                     |
| <b>Uniprot ID</b>          | <b>P62491</b>                                       |
| <b>Dilution Range</b>      | WB: 1:500 - 1:2000IHC: 1:50 - 1:100IF: 1:50 - 1:100 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                     |




Western blot analysis of extracts of var...



Immunofluorescence analysis of C6 cells

...



Immunofluorescence analysis of HeLa cell...