


Product Datasheet

GSK3B Antibody (orb1258847)

Description	GSK3B Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Recombinant protein of human GSK3B.
Target	GSK3B
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IP: 1:50 - 1:100IHC: 1:50 - 1:100
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 47kDa
Uniprot ID	P49841
Dilution Range	WB: 1:500 - 1:2000IP: 1:50 - 1:100IHC: 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...