

Product Datasheet

ADARB1 Antibody (orb1258882)

Description	ADARB1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-380 of human ADARB1 (NP_001103.1).
Target	ADARB1
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 81kDa
Uniprot ID	P78563
Dilution Range	WB: 1:500 - 1:2000
Expiration Date	12 months from date of receipt.