


Product Datasheet


H3K36ac Antibody (orb1259169)

Description	H3K36ac Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	A synthetic acetylated peptide around K36 of human Histone H3 (NP_003484.1).
Target	H3K36ac
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IF: 1:50 - 1:100
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16695
Dilution Range	WB: 1:500 - 1:2000IF: 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of C6 cells

...



Immunofluorescence analysis of NIH/3T3

C...