



Product Datasheet

MAGE-1 antibody (orb125943)

Description	Rabbit Polyclonal to MAGE-1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8WWH6 , P43355
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from H...



Immunohistochemical analysis of formalin...