


Product Datasheet

FN3KRP Antibody (orb1259473)

Description	FN3KRP Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-309 of human FN3KRP (NP_078895.2).
Target	FN3KRP
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:200 - 1:2000
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 35kDa
Uniprot ID	Q9HA64
Dilution Range	WB: 1:200 - 1:2000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...