


Product Datasheet

SLC12A5 Antibody (orb1259490)

Description	SLC12A5 Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 945-1025 of human SLC12A5 (NP_001128243.1).
Target	SLC12A5
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:200 - 1:2000
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 135kDa
Uniprot ID	Q9H2X9
Dilution Range	WB: 1:200 - 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mou...