


Product Datasheet

WARS2 Antibody (orb1259530)

Description	WARS2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 211-360 of human WARS2 (NP_056651.1).
Target	WARS2
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:200 - 1:2000
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 50kDa
Uniprot ID	Q9UGM6
Dilution Range	WB: 1:200 - 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...