



Product Datasheet

CAS9 Antibody (orb1259711)

Description	CAS9 Antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant protein of streptococcus pyogenes CAS9
Target	CAS9
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:5000 - 1:10000IF: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 150kDa
Dilution Range	WB: 1:5000 - 1:10000IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of wil...



Immunofluorescence analysis of HeLa cell...