


Product Datasheet


FGA Antibody (orb1259907)

Description	FGA Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-339 of human FGA (NP_068657.1).
Target	FGA
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	IF: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02671
Dilution Range	IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of C6 cells

...



Immunofluorescence analysis of NIH/3T3

C...



Immunofluorescence analysis of U-2 OS

ce...