



Product Datasheet

MT-ND1 antibody (orb126000)

Description	Rabbit Polyclonal to MT-ND1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P03886
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from A...



Immunohistochemical analysis of formalin...