



Product Datasheet

RBBP4 Antibody (orb1260178)

Description	RBBP4 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, IF, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-425 of human RBBP4 (NP_005601.1).
Target	RBBP4
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IF: 1:50 - 1:200ChIP: 1:20 - 1:50
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 48kDa
Uniprot ID	Q09028
Dilution Range	WB: 1:500 - 1:2000IF: 1:50 - 1:200ChIP: 1:20 - 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Chromatin immunoprecipitation analysis o...