



Product Datasheet

IGF2BP1 Antibody (orb1260345)

Description	IGF2BP1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 480 to the C-terminus of human IGF2BP1 (NP_006537.3).
Target	IGF2BP1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IF: 1:50 - 1:200IP: 1:20 - 1:50
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 70kDa
Uniprot ID	Q9NZI8
Dilution Range	WB: 1:500 - 1:2000IF: 1:50 - 1:200IP: 1:20 - 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunoprecipitation analysis of 200 ug e...