



Product Datasheet

CTCF Antibody (orb1260507)


Description	CTCF Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IHC, IP, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human CTCF (NP_006556.1).
Target	CTCF
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:1000 - 1:2000IHC: 1:50 - 1:200IF: 1:50 - 1:200IP: 1:50 - 1:100ChIP: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 140kDa
Uniprot ID	P49711
Dilution Range	WB: 1:1000 - 1:2000IHC: 1:50 - 1:200IF: 1:50 - 1:200IP: 1:50 - 1:100ChIP: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of 293...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...