


Product Datasheet

ACTL7B Antibody (orb1260567)

Description	ACTL7B Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-380 of human ACTL7B (NP_006677.1).
Target	ACTL7B
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:1000 - 1:2000
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 45kDa
Uniprot ID	Q9Y614
Dilution Range	WB: 1:1000 - 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...