


Product Datasheet

CYP2D6 Antibody (orb1260652)

Description	CYP2D6 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-230 of human CYP2D6 (NP_000097.3).
Target	CYP2D6
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:100 - 1:500 IHC: 1:20 - 1:100
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 60kDa
Uniprot ID	P10635
Dilution Range	WB: 1:100 - 1:500 IHC: 1:20 - 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...