



Product Datasheet

POLR1E Antibody (orb1260783)

Description	POLR1E Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-419 of human POLR1E (NP_071935.1).
Target	POLR1E
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IP: 1:50 - 1:100
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 56kDa
Uniprot ID	Q9GZS1
Dilution Range	WB: 1:500 - 1:2000IP: 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunoprecipitation analysis of 200 ug e...