


Product Datasheet

BRWD1 Antibody (orb1260803)

Description	BRWD1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human BRWD1 (NP_001007247.1).
Target	BRWD1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:1000 - 1:3000
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 263kDa
Uniprot ID	Q9NSI6
Dilution Range	WB: 1:1000 - 1:3000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...