



Product Datasheet

IL7R Antibody (orb1260981)

Description	IL7R Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-239 of human CD127 (NP_002176.2).
Target	IL7R
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IF: 1:20 - 1:50
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 75kDa
Uniprot ID	P16871
Dilution Range	WB: 1:500 - 1:2000IF: 1:20 - 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunofluorescence analysis of HeLa cell...