


Product Datasheet

Bestrophin 4 antibody (orb1261)

Description	Bestrophin 4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Bestrophin 4 (151-250/473 aa)
Target	Bestrophin 4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	53 kDa
Uniprot ID	Q9C0D5
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: HepG2 (Human) Cell; HeLa (Huma...)