


Product Datasheet

PNMAL1 Antibody (orb1261072)

Description	PNMAL1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-439 of human PNMAL1 (NP_060685.2).
Target	PNMAL1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:1000
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 48kDa
Uniprot ID	Q86V59
Dilution Range	WB: 1:500 - 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...