



## Product Datasheet

### TPM3 Antibody (orb1261086)

|                            |                                                                                                                  |
|----------------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | TPM3 Antibody                                                                                                    |
| <b>Species/Host</b>        | Rabbit                                                                                                           |
| <b>Reactivity</b>          | Human, Mouse                                                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                                                                     |
| <b>Tested Applications</b> | IF, WB                                                                                                           |
| <b>Immunogen</b>           | Recombinant fusion protein containing a sequence corresponding to amino acids 1-284 of human TPM3 (NP_689476.2). |
| <b>Target</b>              | TPM3                                                                                                             |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.                                                               |
| <b>Form/Appearance</b>     | Liquid                                                                                                           |
| <b>Concentration</b>       | batch dependent                                                                                                  |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                                                                      |
| <b>Note</b>                | For research use only                                                                                            |
| <b>Application notes</b>   | WB: 1:500 - 1:2000IF: 1:50 - 1:200                                                                               |
| <b>Isotype</b>             | IgG                                                                                                              |
| <b>Clonality</b>           | Polyclonal                                                                                                       |
| <b>MW</b>                  | Observed: 38kDa                                                                                                  |
| <b>Uniprot ID</b>          | <b>P06753</b>                                                                                                    |
| <b>Dilution Range</b>      | WB: 1:500 - 1:2000IF: 1:50 - 1:200                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                  |



Immunofluorescence analysis of HeLa cell...



Western blot analysis of extracts of var...