


## Product Datasheet

### MRPS18B antibody (orb12613)

|                            |                                                           |
|----------------------------|-----------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to MRPS18B                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                         |
| <b>Conjugation</b>         | Unconjugated                                              |
| <b>Tested Applications</b> | IHC, WB                                                   |
| <b>Immunogen</b>           | Recombinant fusion protein of human MRPS18B(NP_054765.1). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.          |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.               |
| <b>Note</b>                | For research use only                                     |
| <b>Isotype</b>             | IgG                                                       |
| <b>Clonality</b>           | Polyclonal                                                |
| <b>MW</b>                  | 35kDa                                                     |
| <b>Uniprot ID</b>          | <b>Q9Y676</b>                                             |
| <b>Dilution Range</b>      | WB: 1:500-1:2000                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                           |



Western blot analysis of extracts of var...