



Product Datasheet

CSRP1 Antibody (orb1261364)

Description	CSRP1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-193 of human CSRP1 (NP_004069.1).
Target	CSRP1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IP: 1:50 - 1:200RIP: 1:20 - 1:50
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 21kDa
Uniprot ID	P21291
Dilution Range	WB: 1:500 - 1:2000IP: 1:50 - 1:200RIP: 1:20 - 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunoprecipitation analysis of 200ug ex...