


Product Datasheet

SMUG1 Antibody (orb1261625)

Description	SMUG1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human SMUG1 (NP_055126.1).
Target	SMUG1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:200 - 1:1000
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 30kDa
Uniprot ID	Q53HV7
Dilution Range	WB: 1:200 - 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...