


Product Datasheet

HRH4 Antibody (orb1261638)

Description	HRH4 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 194-304 of human HRH4 (NP_067637.2).
Target	HRH4
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:1000 - 1:4000
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 44kDa
Uniprot ID	Q9H3N8
Dilution Range	WB: 1:1000 - 1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...