


Product Datasheet

GP1BB Antibody (orb1261666)

Description	GP1BB Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-147 of human GP1BB (NP_000398.1).
Target	GP1BB
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:200 - 1:1000
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 30kDa
Uniprot ID	P13224
Dilution Range	WB: 1:200 - 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...