


Product Datasheet

LRRC59 Antibody (orb1261714)

Description	LRRC59 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-244 of human LRRC59 (NP_060979.2).
Target	LRRC59
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:1000 - 1:4000
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 37kDa
Uniprot ID	Q96AG4
Dilution Range	WB: 1:1000 - 1:4000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...