

Product Datasheet

HIAT1 Antibody (orb1261950)

Description

HIAT1 Antibody

Species/Host

Rabbit

Reactivity

Hamster, Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This HIAT1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 454-483 amino acids from the C-terminal region of human HIAT1.

Target

HIAT1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

53 kDa

Uniprot ID

[Q96MC6](#)

NCBI

[Q96MC6](#)

Dilution Range

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

Expiration Date

12 months from date of receipt.

95

55

36

28

Western blot analysis of HIAT1 Antibody ...

m.CHO

95


72

55

36

28

Western blot analysis of HIAT1 Antibody ...



HIAT1 Antibody IHC analysis in formalin ...