

## Product Datasheet

### DDX6 Antibody (orb1261992)

## Description

DDX6 Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

IHC-P, WB

### Applications

### Immunogen

This DDX6 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 348-375 amino acids from the Central region of human DDX6.

### Target

DDX6

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

54 kDa

### Uniprot ID

[P26196](#)

### NCBI


[P26196](#)

### Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

### Expiration Date

12 months from date of receipt.



Western blot analysis in HeLa, MCF-7 cel...



DDX6 Antibody IHC analysis in formalin f...