

Product Datasheet

MAP4K3 Antibody (orb1262013)

Description

MAP4K3 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This MAP4K3 antibody is generated from rabbits immunized with a recombinant protein of human MAP4K3.

Target

MAP4K3

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

101 kDa

Uniprot ID

[Q8IVH8](#)

NCBI


[Q8IVH8](#)

Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis of MAP4K3 Antibody...



MAP4K3 Antibody IHC analysis in formalin...



Flow cytometric analysis of HepG2 cells

...