

## Product Datasheet

### ATP6V1A Antibody (orb1262056)

## Description

ATP6V1A Antibody

## Species/Host

Rabbit

## Reactivity

Human

## Conjugation

Unconjugated

## Tested

IHC-P, WB

## Applications

## Immunogen

This ATP6V1A antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 441-468 amino acids from the Central region of human ATP6V1A.

## Target

ATP6V1A

## Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

## Form/Appearance

Liquid

## Concentration

batch dependent

## Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

## Note

For research use only

## Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

## Isotype

Rabbit Ig

## Clonality

Polyclonal

## MW

68 kDa

## Uniprot ID

**P38606**

## NCBI



**P38606**

## Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

## Expiration Date

12 months from date of receipt.



Western blot analysis of ATP6V1A Antibod...



ATP6V1A Antibody IHC analysis in formali...