

Product Datasheet

SLC25A37 Antibody (orb1262089)

Description

SLC25A37 Antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications

Immunogen

This SLC25A37 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 309-338 amino acids from the C-terminal region of human SLC25A37.

Target

SLC25A37

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IF starting dilution is: 1:25For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

37 kDa

Uniprot ID

[Q9NYZ2](#)

NCBI



[Q9NYZ2](#)

Dilution Range


For IF starting dilution is: 1:25For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Immunofluorescent analysis of HeLa cells...



Western blot analysis in K562 cell line ...



SLC25A37 Antibody IHC analysis in formalin...