





Product Datasheet


ABHD14B Antibody (orb1262095)

Description	ABHD14B Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	FC, IHC-P, WB	Western blot analysis in A2058, T47D cel...
Immunogen	This ABHEB antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 150-179 amino acids from the C-terminal region of human ABHEB.	
Target	ABHD14B	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	ABHEB Antibody immunohistochemistry anal...
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	22 kDa	
Uniprot ID	Q96IU4	
NCBI	Q96IU4	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	

Western blot analysis in A2058, T47D cel...



ABHEB Antibody immunohistochemistry anal...



Flow cytometric analysis of A2058 cells