

Product Datasheet

HFE2 Antibody (orb1262097)

Description

HFE2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, WB

Applications
Immunogen

This HFE2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 308-338 amino acids from the C-terminal region of human HFE2.

Target

HFE2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


45 kDa

Uniprot ID
Q6ZVN8
NCBI
Q6ZVN8
Dilution Range


For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.



Western blot analysis of HFE2 Antibody i...



Flow cytometric analysis of HepG2 cells

...