

## Product Datasheet

### CHMP4B Antibody (orb1262174)

## Description

CHMP4B Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

FC, IHC-P, WB

### Applications

### Immunogen

This CHMP4B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 51-79 amino acids from the N-terminal region of human CHMP4B.

### Target

CHMP4B

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

25 kDa

### Uniprot ID

[Q9H444](#)

### NCBI

[Q9H444](#)


### Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50


### Expiration Date

12 months from date of receipt.


Sk-Br-3  
95



Western blot analysis of CHMP4B Antibody...



Formalin-fixed and paraffin-embedded hum...



Flow cytometric analysis of MCF-7 cells

...