

Product Datasheet

IL2RG Antibody (orb1262210)

Description

IL2RG Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, WB

Applications

Immunogen

This IL2RG antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 76-101 amino acids from the N-terminal region of human IL2RG.

Target

IL2RG

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IF starting dilution is: 1:25For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

42 kDa

Uniprot ID

[P31785](#)

NCBI



[P31785](#)

Dilution Range


For IF starting dilution is: 1:25For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Immunofluorescent analysis of 4% parafor...



Western Blot at 1:2000 dilution Lane 1: ...



Western blot analysis of lysate from Jur...