

Product Datasheet

FBXL5 Antibody (orb1262231)

Description
FBXL5 Antibody
Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This FBXL5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 86-115 amino acids from the N-terminal region of human FBXL5.

Target

FBXL5

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


79 kDa

Uniprot ID
Q9UKA1
NCBI
Q9UKA1
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis of FBXL5 Antibody ...



Western blot analysis of FBXL5 Antibody ...



FBXL5 Antibody IHC analysis in formalin ...