

Product Datasheet

TRERF1 Antibody (orb1262253)

Description

TRERF1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, WB

Applications

Immunogen

This TREF1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1136-1164 amino acids from the C-terminal region of human TREF1.

Target

TRERF1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IF starting dilution is: 1:10~50For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

132 kDa

Uniprot ID

Q96PN7

NCBI



Q96PN7

Dilution Range


For IF starting dilution is: 1:10~50For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Fluorescent image of MCF-7 cell stained ...



Western blot analysis of TREF1 Antibody ...



Flow cytometric analysis of K562 cells (...)