

Product Datasheet

HSD17B12 Antibody (orb1262256)

Description

HSD17B12 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This HSD17B12 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 126-155 amino acids from the Central region of human HSD17B12.

Target

HSD17B12

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

34 kDa

Uniprot ID

Q53GQ0

NCBI


Q53GQ0

Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis in A431, SiHa cell...



Formalin-fixed and paraffin-embedded hum...



Flow cytometric analysis of A2058 cells

...