

## Product Datasheet

### EXTL2 Antibody (orb1262309)

**Description**

EXTL2 Antibody

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse

**Conjugation**

Unconjugated

**Tested**

WB

**Applications**
**Immunogen**

This EXTL2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 269-295 amino acids from the C-terminal region of human EXTL2.

**Target**

EXTL2

**Preservatives**

Supplied in PBS with 0.09% (W/V) sodium azide.

**Form/Appearance**

Liquid

**Concentration**

batch dependent

**Storage**

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

For WB starting dilution is: 1:1000

**Isotype**

Rabbit Ig

**Clonality**

Polyclonal

**MW**

37 kDa


**Uniprot ID**
**Q9UBQ6**
**NCBI**
**Q9UBQ6**
**Dilution Range**

For WB starting dilution is: 1:1000


**Expiration Date**

12 months from date of receipt.

MDA-MB231



Western blot analysis of EXTL2 Antibody ...



Western blot analysis of EXTL2 Antibody ...