

Product Datasheet

TSPAN6 Antibody (orb1262312)

Description

TSPAN6 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

This TSPAN6 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 176~205 amino acids from the C-terminal region of human TSPAN6.

Target

TSPAN6

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

28 kDa

Uniprot ID

[O43657](#)

NCBI


[O43657](#)

Dilution Range


For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis of TSPAN6 Antibody...



Flow cytometry analysis of HeLa cells (b...



Western blot analysis of TSPAN6 using ra...